
BioBundle, complete cultivation systems

A BioBundle is a complete bioreactor system, equipped with all necessary components and is ready to use "out-of-the-box". No detail is overlooked: the system is complete with silicone tubing, sample bottles and a "starter kit" including spare parts. The BioBundle is easy to set-up, requires no special skills or tools, easy to learn and easy to operate. BioXpert Lite Software for data acquisition is included. Select one or more of the optional add-on packs to customize your BioBundle. The BioBundle provides a unique



combination of ease of use and sophistication, reducing the time to start-up a process. The system is equipped with the intelligent and powerful process controller: easy and intuitive to operate, such that the user manual might not be needed.

The process controller has control loops for pH, Temperature, Dissolved Oxygen, Foam/Level and agitation and can apply a combination of digital and analog outputs for process control. Actuators such as rotameters, solenoid valves, mass flow controllers, pumps, thermocirculator and others can be controlled. The operator can set P-I-D values, dead-band for pH control, cascade control strategies, and dose monitors for liquid additions. The advanced auto-tuning adaptive control system is part of every bundle and takes the guessing out of PID controller setting. The system can automatically and continuously calculate the best controller settings for every process.



Modularity and configurability

Although the BioBundle has a pre-set configuration, the whole bioreactor system remains modular and configurable. With simple changes in configuration in the BioBundle system can be modified for different applications in a cost-effective way! This applies to changes in process control strategies, in using extra gas or liquids, but also in changing the system into a microbial set-up.

Compact design

The BioBundle is designed to occupy as small a footprint as possible. All accessories such as pumps and gas flow control valves are conveniently integrated in a compact console. The compact design reduces the need for expensive lab space.

Software for data acquisition and supervisory control

The BioBundle control system includes an Ethernet connection port to connect to a PC with software for data acquisition or SCADA such as the Applikon BioXpert packages.

BioXpert Lite (data acquisition) is included in this BioBundle. Multiple bioreactors can be connected to one PC for data acquisition.

Optional BioXpert packages are BioXpert 2 and BioXpert W7: SCADA of multiple bioreactors, including cGMP production applications.

Cost Saving

The BioBundle is a pre-packed system. The assembling of the BioBundles in series provides a cost saving which is for the benefit of the customer.

Summarising

Over all some of the benefits Applikon bioreactors provide include:

- dependable and reliable operation = greater productivity and yield
- modular design = less cost in the future
 - (being able to use same equipment for different applications)
- less downtime due to maintenance = greater economy in operation, higher yield to cost ratio
- simple operation = less operator training required, the product

Specifications

	MiniBioBundle	MiniBioBundle	BioBundle	BioBundle	Appliflex
	Cell Culture	Microbial	Cell Culture	Microbial	BioBundle
Control system	my-Control	my-Control	ez-Control	ez-Control	ez-Control
Total Volume	250ml, 500 ml,	250ml, 500 ml,	1L, 2L, 3L, 5L,	1L, 2L, 3L, 5L,	10L, 20L, 50L
	1000ml	1000ml	7L, 15L, 20L	7L, 15L, 20L	
Working Volume	200ml, 400ml,	200ml, 400ml,	0.9L, 1.7L, 2.7L,	0.9L, 1.7L, 2.7L,	5L, 10L, 25L
	800ml	800ml	3.2L, 5.4L, 12L, 16L	3.2L, 5.4L, 12L, 16L	
Agitator	Lipsealed with	Lipsealed with	Lipsealed with	Lipsealed with	Mixing by
	marine impeller	Rushton impellers	marine impeller	Rushton impellers	rocking motion
Aeration	Air supply via	Air supply via	Air and Oxygen	Air and Oxygen	Air and Oxygen
	sparger	sparger	supply via sparger	supply via sparger	supply via overlay
				and overlay	
Exhaust gas	Optional gas outlet	Gas outlet	Gas outlet	Gas outlet	
		condenser	condenser	condenser	condenser
Sampling	Sample pipe	Sample pipe	Sample pipe	Sample pipe	Sample line
	included	included	included	included	included
	Sample system				
	optional	optional	optional	optional	optional
pH	Measurement via				
	pH sensor control				
	via liquid alkali	via liquid alkali or	via liquid alkali	via liquid alkali and	via liquid alkali
	pump and CO ₂	acid addition pump	pump and CO ₂	acid addition pump	pump and CO ₂
	gas supply		gas supply		gas supply
Temperature	Measurement via				
	Pt-100	Pt-100	Pt-100	Pt-100	Pt-100
	Heating only via	Heating and	Heating only via	Heating via heating	Heating only via
	heating blanket	cooling via	heating blanket	blanket cooling by	heating blanket
	<u> </u>	Peltier system	Ü	cold water in heat	G
		•		exchanger	
Dissolved Oxygen	Measurement via				
	DO ₂ sensor control	DO ₂ sensor contro			
	via Air and O ₂ gas	via Air supply and	via Air and O ₂ gas	via Air and O ₂	via Air and O ₂ gas
	supply	agitation speed	supply	supply and agitation	supply
				speed	
Foam	Option	Measurement via	Option	Measurement via	Option
		Foam sensor		Foam sensor	
		control via anti-foam	1	control via anti-foam	
		addition pump		addition pump	
Level	Option	Option	Option		Option
Liquid additions	4 fixed in topplate	4 fixed in topplate	3 ports in triple	3 ports in triple inlet,	One inoculum and
	and 1 septum port	and 1 septum port	inlet and one	one inoculum port	alkali inlet.
		•	inoculum port.	1 and septum port.	
	One liquid storage	Two liquid storage	One liquid storage	Three liquid storage	One liquid storage
	system included	systems included	system included	systems included	system included
Start-up kit	Included	Included	Included	Included	Included
BioXpert Lite software	Included	Included	Included	Included	Included
	-				